

Low sulphur in various samples
Sensitivity < 20 ppb
High stability, low maintenance
Factory tuned

APPLICATION NOTE 212WA1012E

Low sulphur analyser

ASTM D6228, D5303 D5504,
D3328, D4735, D5623,
D7011

UOP 791
ISO 19739

G-A-S offers custom configured GC analysers for complex separations, data processing and reporting. We have over 35 years of experience in designing and building turnkey analysers for many application fields. Our analysers are designed to meet many accepted standard methods (like GPA, ASTM, UOP, ISO, etc.) in the Oil and Gas industry. The efficient configurations are based on proven GC technology, resulting in robust instruments with an optimal return on investment.

The determination of sulphur components in natural gas and various hydrocarbon streams is of utmost importance for several reasons. Sulphur is known to be very harmful to expensive catalysts involved in downstream hydrocarbon processing. Sulphur components are corrosive to equipment, and sulphur amount in finished fuel products is lowered worldwide to reduce emission. GC-(P)FPD -Gas Chromatography with (Pulsed) Flame Photometric Detector- has proven to be a very stable and reliable solution for sulphur determination in many laboratories.



Figure 1: Low sulphur analyser with optional Sample Securitiser

Gaseous samples like natural gas and ethylene are injected using GSV (Gas Sampling Valve), according to ASTM D6228, D5303, D5504, UOP791 and ISO 19739.

LPG samples are analysed by GSV after vaporising or by LSV (Liquid Sampling Valve) using direct liquid injection. The optional Sample Securitiser is used for guaranteed quantitative injection by LSV. See figure 1 and our dedicated application note for more information.

Liquid samples like finished fuels are injected using SSL (Split injector) or PTV (Programmable Temperature Vaporiser), according to ASTM D5623.

Sulphur components very easy adsorb on active surfaces, therefore Sulfinert deactivation is applied to the full sample path, which is vital for accurate quantification.




Figure 2: Low sulphur analysis in propylene (LSV-Rt SilicaBond-PFPD; 1 ppm per component)




Figure 3: Low sulphur in natural gas (GSV-Rtx 1,FPD)




Figure 4: Low sulphur in diesel fuel (SSL-Rtx 1-PFPD)




Figure 5. Calibration standard of sulphur components in N₂ (5 ppm each) (GSV-Rtx 1-PFPD)




Figure 6. Linearity curve of H₂S: $r^2 = 0.9997$ (100ppb - 5 ppm) ($r^2 > 0.995$ for all components; PFPD, lineariser mode)

GC-PFPD ID	H ₂ S Area	COS Area	CH ₃ SH Area	C ₂ H ₅ SH Area	i-C ₃ H ₇ SH Area	C ₃ H ₇ SH Area	THI Area	THTHI Area
PFPD_071	7653680	7718514	6026713	6764998	7612207	7096988	9317300	9228699
PFPD_072	7575351	7695225	5928568	6794046	7631111	7204149	9296086	9213385
PFPD_073	7648397	7698897	6013745	6725392	7641811	7193838	9352864	9246370
PFPD_074	7608142	7713938	6034135	6695520	7598043	7161945	9327279	9249580
Min:	7575351	7695225	5928568	6695520	7598043	7096988	9296086	9213385
Max:	7653680	7718514	6034135	6794046	7641811	7204149	9352864	9249580
Mean:	7621393	7706144	6000790	6744989	7620793	7164230	9323382	9234509
Std Dev:	36820	10646	48880	43352	19489	48294	23569	16811
%RSD:	0.48	0.14	0.81	0.64	0.26	0.67	0.25	0.18

Figure 7. Repeatability at 5 ppm concentration level (GSV-PFPD).

Component	LOD (ppb)
H ₂ S	17
COS	18
CH ₃ SH	25
C ₂ H ₅ SH	27
i-C ₃ H ₇ SH	39
C ₃ H ₇ SH	30
THI	22
THTHI	30

Figure 8. LOD (s/n = 3 * noise) (PFPD)

Component	LOD (ppb)	
	PFPD	FPD
H ₂ S	33	176
COS	31	158
CS ₂	35	196
C ₂ H ₆ S	41	206
CH ₄ S	18	90

Figure 9. Comparison of PFPD and FPD sensitivity. Conditions: column: Rtx-1, 60 m * 0.32 mm ID, 5 µm; sample loop: 250 µl; column flow: 2 ml/min; splitflow: 5 ml/min; test standard: 500ppb per component; (LOD based on s/n=3)



Figure 10. Optional calibrator module for generating multi level calibration standards, using permeation tubes or using dilution of a calibration gas.

Very sensitive. High selectivity

PFPD provides high sensitivity with LODs below 100ppb. Even lower limits (< 20ppb) are reached when sulphur components are well-separated from the hydrocarbon matrix. Good selectivity is offered as well; hydrocarbon components show no positive response. However co-elution with (very) high amounts of hydrocarbons decreases sensitivity and therefore chromatographic separation is optimised for each standardised method.

PFPD or FPD ?

PFPD and FPD are both available; the first one shows best sensitivity (20-30ppb) and selectivity while the latter provides a more cost effective solution, with detection limits of 100-200 ppb. See figure 8 and 9; note that LODs for PFPD are different in both figures due to the used conditions.

High stability. Easy operation compared to SCD

There is no need for frequent calibration and maintenance thanks to a high long term stability. (P)FPD is easier to operate, is less expensive and the response is more stable compared to the Chemiluminescence detector.

Combining methods

Low sulphur analysis can be combined efficiently with for example refinery gas or natural gas analysis, using a single instrument (figure 11).

Specifications

Standardised methods:

Standard Designation (Committee)	D3328 (D19)	D4735 (D16)	D5303 (D02)	D5504 (D03)	D5623 (D02)	D6228 (D03)	D7011 (D16)
Title	Comparison of Waterborne Petroleum Oils by Gas Chromatography	Determination of Trace Thiophene in Refined Benzene by Gas Chromatography	Trace Carbonyl Sulfide in Propylene by Gas Chromatography	Determination of Sulfur Compounds in Natural Gas and Gaseous Fuels by Gas Chromatography and Chemiluminescence	Sulfur Compounds in Light Petroleum Liquids by Gas Chromatography and Sulfur Selective Detection	Determination of Sulfur Compounds in Natural Gas and Gaseous Fuels by Gas Chromatography and Flame Photometric Detection	Determination of Trace Thiophene in Refined Benzene by Gas Chromatography and Sulfur Selective Detection
Matrix	Petroleum oils recovered from water or beaches compared to oils from suspect sources.	Refined benzene	Propylene	High methane-content gaseous fuels such as natural gas	Light petroleum liquids such as distillates, gasoline motor fuels, and other petroleum liquids with final boiling point of 230 °C or lower.	Gaseous fuels	Refined benzene
Analytes	Hydrocarbons (FID) and associated sulfur compounds (FPD)	Thiophene	Carbonyl sulfide (COS)	Speciated volatile sulfur-containing compounds such as H ₂ S, COS, SO ₂ , thiols, sulfides, thiophenes, etc.	Volatile sulfur-containing compounds; both unidentified and known individual compounds are determined.	Volatile sulfur-containing compounds; both unidentified and known individual compounds are determined.	Thiophene
Range	Qualitative only; no quantitative range is cited	0.5 to 5.0 mg/kg	0.5 to 4.0 mg/kg	0.01 to 1,000 mg/m ³	0.1 to 100 mg/kg (ppm)	0.02 to 20 mg/m ³ (0.014 to 14 ppmv)	0.02 to 2.0 mg/kg
Detector Cited	FID alone, or FID with flame photometric detection	Any flame photometric detector (FPD or PFPD)	A flame photometric detector	FID and SCD in series, or other sulfur specific detector	A sulfur selective detector that produces a linear and equimolar response to sulfur compounds	Any flame photometric detector (FPD or PFPD)	Pulsed flame photometric detector, sulfur chemiluminescence detector, atomic emission detector, or any other sulfur selective detector.
Detector Specification	<ul style="list-style-type: none"> A flame photometric detector is cited when independent, selective detection of HC and S is desired to improve identification of oil Linearization of the sulfur signal is recommended 	Signal-to-noise of at least 2-to-1 for 0.5 mg/kg (ppm) thiophene in benzene (4-µL injection with FPD or 1-µL injection with PFPD)	Signal-to-noise of at least 2-to-1 for 0.1 mg/kg COS in propylene	Sulfur specific detectors other than the SCD can be used provided that they have sufficient sensitivity, respond to all eluted sulfur compounds, do not suffer from interferences, and satisfy quality assurance criteria.	<ul style="list-style-type: none"> Linearity of 10⁴ 5-pg sulfur/second minimum detectability Approximate equimolar response to sulfur No interference from co-eluting hydrocarbons 	Any flame photometric detector calibrated in the sulfur-specific mode is used for this test method	<ul style="list-style-type: none"> Linearity of 10⁷ Minimum detectable level of less than 0.02 mg/kg thiophene in benzene Selectivity of sulfur to carbon greater than 10³ Absence of quenching that affect results
PFPD Suitability	A PFPD with increased sensitivity, equimolar sulfur response, and linear signal provides maximum information about the number, distribution, and relative concentrations of the sulfur species present for positive identification of the suspect oil.	The PFPD provides additional sensitivity, down to 0.02 mg/kg (ppm), when lower concentrations of thiophene in benzene must be measured.	Linear, equimolar response of the PFPD eliminates log/log plots, simplifies calibrations, and enables quantitation of any unknown sulfur peaks that may be present along with the COS.	The PFPD exceeds requirements for "other sulfur specific detectors", including all QA criteria. The two-signal output of the PFPD provides a simultaneous fingerprint of the hydrocarbon matrix along with the sulfur chromatogram.	Linear, equimolar response of the PFPD allows accurate quantitation of individual known and unknown sulfur compounds down to 25-ppb sulfur or lower.	Increased sensitivity of the PFPD allows quantitation of sulfur at concentrations 5 to 10 times lower than an FPD. Linear, equimolar response of the PFPD allows accurate quantitation of the sulfur content in unknowns.	Long-term stability and ease of use make the PFPD well suited for use in production facilities with limited staff, and where low detection limits are required.

- UOP 791 and ISO 19739 are covered as well

- Configuration:** 1-4 channel instrument based on Thermo Trace 1300 GC or CompactGC^{4.0}.
- Detectors:** Pulsed Flame Photometric Detector or Flame Photometric Detector
- Optional:**
- Vaporiser or Sample Securitiser for LPG samples
 - hydrocarbon signal output on PFPD
 - simultaneous detection on FID using effluent splitter
- Sample tubing:** Sulfinert[®] tubing for inert sample path
- Application:** Custom configured analyser for the analysis of sulphur components in various gases, liquefied gases and liquid streams
- Range:** 20 ppb - 100 ppm (PFPD); 200ppb-100ppm (FPD)
- Detection limit**
- <100 ppb; <20 ppb in case of well-separated sulphur components (PFPD)
 - < 200 ppb (FPD)
- Repeatability** < 3 % RSD or better
- Sample requirements:** See our pre-installation guide for additional requirements
- (P)FPD compared to SCD:** (P)FPD is easier to operate and offers better long term stability



Figure 11. Refinery Gas and Low sulphur analysis offered in a single instrument.



Figure 12. PFPD is also available on 19" CompactGC^{4.0} for fast ppb sulphur analysis in several matrices.

GLOBAL™ ANALYSER SOLUTIONS

G·A·S

GAS is an INTERSCIENCE company

GAS is the preferred solution partner of

Thermo SCIENTIFIC

P.O. Box 2148, 4800 CC Breda, The Netherlands

T: +31 (0)76 5411800 F: +31 (0)842 206757 info@gas-site.com www.gas-site.com