



TOGA

Transformer Oil Gas Analyser

- Fully automated Transformer Oil Gas Analyser
- According to ASTM D3612c
- High uptime due to modular injector/detector technology
- High simplicity and convenience by integration with Chromeleon CDS

AN 211WA0620F

GAS offers custom configured GC analysers for many application fields since 40 years. GAS analysers are designed to meet many standardised methods from GPA, ASTM, UOP, ISO, EN and others. The efficient configurations are based on proven GC technology, resulting in robust instruments with an optimal return on investment. Transformer oil is a highly refined mineral oil used in electrical transformers. It has excellent insulating properties, suppresses corona and arcing, and serves as a coolant. In case of electrical errors, the oil breaks down to gases, which identity and content can be related to the type and severity of the electrical fault. ASTM D3612c describes extraction and determination of gases in transformer oil based on head space sampling.



Figure 1. Transformer oil gas analyser with headspace sampling



Figure 2. Schematic diagram TOGA

Transformer Oil Gas Analyser - Principle

An aliquot from the vapor phase above a transformer oil sample in a sealed vial is injected using a headspace autosampler. After H₂, O₂, N₂, CH₄ and CO have passed the first separation column (Rt-UBond), these components are isolated on the second column (Molsieve) by switching valve E1. Next CO₂, C₂H₂, C₂H₄ and C₂H₆ elute from the first column via the needle valve (NV) to the detectors. After valve E1 is switched back to the initial position, H₂, O₂, N₂, CH₄ and CO are detected.

When C_3 and C_4 components need to be analysed as well, a valve and separation column is added to the system. CO and CO_2 are converted to CH_4 using a methaniser, enabling sensitive FID detection at ppm level. (See figure 2).

Diaphragm valves and InstantConnect detector technology result in a robust and flexible instrument.

Headspace sampling

The oil samples are automatically injected using Thermo TriPlus 500 headspace autosampler with 12, 120 or 240 sample position tray. This system uses the closed loop injection principle, so loss of components or false results from ambient gases are avoided. Each sample is individual equilibrated at 70 °C. The required equilibration time is minimised by using the mixer function of the instrument. The headspace sample is injected by (1) pressurisation of the vial (2) expansion to a loop (3) transferring the loop content to the separation columns (see figure 3). To guarantee the expected sample integrity at every analysis, each vial is automatically checked for possible leaks just prior to the loop filling stage. Chromeleon chromatography data system is seamlessly integrated with TriPlus 500 HS autosampler, offering a single point of control for optimal

simplicity and convenience.



Figure 3. TriPlus 500 headspace autosampler; vial pressurising, loop filling and injection

Results

Figure 4 and 5 show the results of 10 and 100 ppm dissolved gases in oil standards on FID and TCD. Figure 6 shows achieved limit of detection per component and applied detector.



Figure 4. FID, 10 ppm Morgan Schaffer oil standard Figure 5. TCD, 100 ppm Morgan Schaffer oil standard Figure 6. LODs TOGA analyser

Specification

Standard methods:	ASTM D3612 part c			
Configuration:	One channel instrument based on Thermo Trace GC1600. using microTCD, methaniser and FID. Automated injection, closed loop principle, using Thermo Triplus 500 (12, 120 or 240 vials of 10, 20 or 22 ml) Carrier gas: Argon			
Application:	Custom configured analyser for the analysis of dissolved gases in transformer oil. Components: H_2 , O_2 , N_2 , CH_4 , CO , CO_2 , C_2H_2 , C_2H_4 , C_2H_6 . C_3 and C_4 components optional.			
Sample requirements:	The oil sample must be offered to the analyser using the appropriate 10 or 20/22 sample vials. Vials are purged with Argon before sampling (see optional sample preparation device (figure 7).			
Analysis Time:	20 minutes.			
Minimum Detectability:	Component	Detection limit GAS (ppm)]	-
	CO ₂	< 1	1	
	Acetylene	< 1	-	
	Ethvlene	< 1	-	
	Ethane	< 1		
	Hydrogen	< 5		
	Oxygen	< 100		
	Nitrogen	< 100	1	
	Methane	< 1		
	СО	< 1]	
Dynamic Range:	4 decades for TCD, 7 decades for FID.7 Figure 7. Optional DGA Revolver Table system for TOGA sample preparation			
Repeatability:	Better than 5% RSD at 100ppm concentration level for all analytes specified, measured over at least 10 consecutive			
	runs.			
Optional:	Additional gas sampling valve for injection of gas samples without use of the autosampler Additonal column and switching valve for C_3/C_4 analysis DGA Revolver Table system for TOGA sample preparation (figure 7; avoiding air interference)			
Data systems:	Chromeleon CI	DS		



GAS is an Interscience company

WWW.GASSITE.COM