

Hydrogen purity analyser

- Green hydrogen testing
- ISO 21087, ISO 14687, EN 17124
- PEM Fuel cell applications and risk assessments

AN31WA0921B

GAS offers custom configured GC analysers for many application fields for over 45 years. GAS analysers are designed to meet many standardised methods from GPA, ASTM, UOP, ISO, EN and others. The efficient configurations are based on proven GC technology, resulting in robust instruments with an optimal return on investment.

Fuel cell electric vehicles are considered to be a very promising alternative to conventional automobiles. They provide zero emissions (only H₂O is produced), high efficiency, long range and short refuelling time. Fuel cell operation relies mainly on platinum catalysts and PEM (proton exchange membrane), both delicate components which are highly sensitive to the quality of the hydrogen fuel.

Various GC & GC-MS methods

International standards ISO 14687 & 21087 address the detailed analysis of impurities in hydrogen. Both standards describe various GC and GC-MS methods for analysing ppb and ppm levels of a range of components. ISO 14687 specifies a restricted list of suggested techniques, while 21087 allows other techniques as well, as long as they are validated. Figure 1 shows an overview of components, maximum allowable levels and GC detectors.

ISO 14687

ISO 14687 describes the following detectors: FID (Flame Ionisation Detector)+methaniser, TCD (Thermal Conductivity Detector), PDD (Pulsed Discharge Detector), SCD (Sulphur Chemiluminescence Detector) and MS (Mass Spectrometer). Components of interest can be pre-concentrated using thermal desorption technique (TD) to obtain lower detection limits. See figure 2, 3 and 4.

Impurity	Limit 14687 21087 (umol/mol)	Spec. 17124 (umol/mol)	ISO 14687 GC method	ISO 21087 GC method (GAS)
Total non-hydrogen gases	300	300	TCD, PDD	PDD
Total hydrocarbons	2 (incl. CH ₄)	2 (excl. CH ₄)	FID, TD-GCMS	MS*, PDD, FID
Methane		100	FID	PDD, FID
Oxygen	5	5	TCD	PDD
Helium	300	300	TCD	TCD
Nitrogen	100	300	TCD, PDD	PDD
Argon	100	300	TCD, PDD	PDD
Carbon dioxide	2	2	PDD or methaniser-FID	PDD or methaniser-FID
Carbon monoxide	0.2	0.2	PDD or methaniser-FID	PDD or methaniser-FID
Total sulphur components	0.004	0.004	SCD TD-SCD	MS*
Formaldehyde	0.01	0.2	FID, PDD, TD-GCMS	MS*
Formic acid	0.2	0.2	TD-GCMS	MS*
Ammonia	0.1	0.1	-	MS* or PDD
Halogenated components	0.05	0.05	ECD, TD-GCMS	MS*



Figure 2. Hydrogen purity analyser with dual GC oven, FID, TCD, PDD, mass spectrometer and thermal desorption option

Figure 1. Impurities, maximum allowable levels and suggested GC detectors according to ISO 14687 & ISO 21087. EN 17124 fuel spec. MS*: ISQ 7000 with AEI (Advanced Electron Ionisation)




Figure 4. Chromatogram SCD. 50 ppm calibration standard.



Figure 3. Low level sulphur analyser with SCD. Limit of detection SCD is < 10 ppb. A pre-concentration device (TD) can be added to increase sensitivity.

ISO 21087

Figure 5 (see next page) shows the hydrogen purity analyser based on ISO 21087, suggested by GAS. Compared to ISO 14687, it offers (1) lower detection limits, (2) a larger number of analysed components and (3) a less complex analyser configuration. PDD and TCD measure permanent gases, while other components, including ammonia, sulphur components, formic acid, formaldehyde and halogenated compounds, are analysed by MS. The AEI (Advanced Electron Ionisation) source provides detection limits below 0.001 vppm without pre-concentration (TD). Ammonia can be analysed by either MS or PDD, see figure 6, 7, 8, 9 and 10 for chromatograms.

EN 17124

EN 17124 specifies quality characteristics of hydrogen fuel and the corresponding quality assurance (see figure 1).




Figure 9. Ammonia using Pulsed Discharge Detector. 13.8 ppm calibration standard in N₂. Limit of detection < 100 ppb.




Figure 6. Permanent gases using PDD (10 ppm calibration standard)




Figure 10. Ammonia using mass spectrometer with Advanced Electron Ionisation option




Figure 7. Halogenated components using mass spectrometer with Advanced Electron Ionisation option




Figure 8. H₂S and COS. MS-AEI detection. 2.2 - 56.6 ppb concentration levels (without pre-concentration)

Specification

Application:	Hydrogen purity analyser, fuel cell PEM.
Standardised methods:	ISO 14687, ISO 21087, EN 17124
GC instrument:	ISO 14687: GC 1300 with various detectors like TCD, FID (+ methaniser), PDD, SCD, PFPD, MS, Thermal Desorption
	ISO 21087: GC 1300 with MS-AEI (Advanced Electron Ionisation), PDD and TCD
Options:	Stop flow valve (Back) pressure regulator Heated sample interface Selector valves (optionally heated), 4-16 sample inlets or multiple, sequence programmed by Chromeleon datasystem. Sample pump Partical filters
Analysed components:	Permanent gases, hydrocarbon, sulphur components, halogenated components, formaldehyde, formic acid, ammonia
Minimum detectability:	Dependant on GC detector: PDD < 50 ppb MS-AEI < 1 ppb TCD < 5 ppm SCD < 10 ppb
Analysis time:	Dependant on used methods
Sample requirements:	See our pre-installation guide



Figure 5. Impurities in hydrogen analyser using PDD, TCD and MS-AEI



Turnkey customised
GC & GC-MS solutions



Xpert center for chromatography

Expert & education centre
Learn from the Xperts!



Fully automated solutions
for sample preparation

GAS, IS-X, IS-X Academy & SampleQ are INTERSCIENCE brandings