



Carbon Capture, Utilisation & Storage gas analyser

- For assuring CO₂ purity
- Combined GC and FTIR methodology
- 19" rack size; small footprint,
- Fast results (< 3 minutes)
- Highly sensitive (LOD ppm-ppb range)

Get ready for **tomorrow's analytics**

CCUS (Carbon Capture, Utilisation and Storage) is considered a critical technology for achieving significant reductions in CO₂ emissions, particularly in sectors that are difficult to decarbonise, such as heavy industry and fossil fuel-based power generation. It is also seen as an important component of strategies to achieve net-zero emissions and mitigate the impacts of climate change. Understanding and managing impurities in CO₂ is crucial for the efficiency, safety, cost-effectiveness, and environmental compliance of CCUS projects.

Combined GC and FTIR methodology

To ensure specification, CO₂ must be checked for several possible impurities. Gas Chromatography (GC) is known for its specific detection capabilities using dedicated separation columns tuned to the required components of interest. Additionally, a wide range of detectors is available, offering either a universal response to several components or selective detection of specific component groups, thereby isolating them from the matrix. With CompactGC^{4.0} (figure 2) fast run times are provided, usually below 3 minutes.

FTIR (Fourier Transform Infrared Spectroscopy) complements GC by providing detailed molecular identification through the detection of specific functional groups and molecular vibrations. It offers rapid, precise analysis of chemical structures, enhancing GC's separation capabilities with additional molecular information using spectral data.

Instrument configurations

CompactGC^{4.0} is equipped with 4 analytical channels, using the following detectors:

- Thermal Conductivity Detector (TCD) for detecting permanent gases and light hydrocarbons.
- Flame Ionisation detectors (2) for analysing hydrocarbons and other carbon containing species.
- PFPD for selective and sensitive detection of sulphur components.

(Configuration may vary due to different CCUS specification)

Fourier Transform Infrared (FTIR) Spectroscopy works by measuring how a sample absorbs infrared (IR) light at different wavelengths. A broad-spectrum IR source passes through the sample, and the transmitted light is collected. A interferometer modulates the IR light, creating an interferogram, which is then converted into an absorption spectrum using Fourier Transform. The resulting spectrum reveals characteristic absorption peaks corresponding to the molecular vibrations of chemical bonds in the sample (figure 3).

FTIR is highly effective for detecting gases like water vapour, ammonia, carbon monoxide, carbon dioxide, formaldehyde and others because of their unique IR absorption fingerprint due to molecular vibrations. FTIR can detect multiple gases simultaneously with high sensitivity. It is widely used in environmental monitoring, industrial safety, and atmospheric studies for gas analysis.




Figure 1 Diagram CCUS gas analyser



Figure 2 CompactGC^{4.0}




Figure 3 FTIR principle

Results

Figures 4, 5, and 6 show example chromatograms of the CompactGC4.0 using TCD, FID, and PFPD detectors, respectively. Figures 7, and 8 display calibration of the FTIR gas analyser. Figures 9 and 10 present the analysed components and their limits of detection.




Figure 4 Permanent gas on TCD




Figure 5 Sulphur components on PFPD




Figure 6 Hydrocarbons on FID




Figure 7 FTIR calibration for Acetaldehyde




Figure 8
FTIR linearity for NO

Specifications

Configuration:

- 4 channel CompactGC^{4.0} equipped with:
 - TCD (Thermal Conductivity Detector)
 - FID (Flame Ionisation detector) (2)
 - PFPD (Pulsed Flame Photometric Detector)
- FTIR gas analyser equipped with DTGS detector

Application:

Carbon Capture, Utilisation and Storage gas analyser; analysis of impurities in CO₂

Analysed components:

- CompactGC: H₂, Ar, O₂, CH₄, CO, hydrocarbons, sulphur components.
- FTIR: CO₂, CO, CH₄, H₂O, methanol, ethanol, acetaldehyde, ammonia, formaldehyde, NO, NO₂

Optional:

- Pressure reducer; back-pressure regulator
- Sample pump
- Stop flow valve
- Stream selector

Sample tubing:

Sulfinert[®] tubing for inert sample path (sulphur and other polar components)

Sample requirements:

See our pre-installation guide for additional requirements

Analysis time:

3 minutes, depending on analysed components

Minimum detectability:

See figure 9 & 10

Data systems:

Chromeleon

FTIR		
Gas	unit	LOD (3σ)
Acetaldehyde	ppm	0.002
Ammonia	ppm	0.01
Carbon dioxide	%	n/a
Carbon monoxide	ppm	0.07
Ethanol	ppm	0.05
Formaldehyde	ppm	0.09
Methane	ppm	n/a
Methanol	ppm	0.07
Nitric oxide	ppm	0.05
Nitrogen dioxide	ppm	0.003
Nitrous oxide	ppm	0.10
Water ppm	ppm	0.41

Figure 9 Analysed components and limits of detection (FTIR).

CompactGC			
Gas	unit	LOD (3*noise)	Detector
Hydrogen	%	< 0.1	TCD
Argon	ppm	< 10	TCD
Oxygen	ppm	< 10	TCD
Nitrogen	ppm	< 10	TCD
Methane	ppm	< 10	TCD
Hydrocarbons *	ppm	< 1 *	FID
Sulphur components	ppm	< 0.1	PFPD
Amines *	ppm	< 1 *	FID
* dependant on component			

Figure 10 Analysed components and limits of detection (CompactGC^{4.0})

powered by **interscience**



Turnkey customised
GC & GC/MS solutions



Expert & education centre
Learn from the Xperts!



Fully automated solutions
for sample preparation

GAS, IS-X & SampleQ are Interscience brands

Find out more on www.gassite.com