

# **High Purity Analyser**

- ASTM D8098
- ppb permanent gases in
   C2 and C3 hydrocarbon streams
- Using Pulsed Discharge Detector (PPD, PDHID)

Get ready for tomorrow's analytics



# Trace analysis of permanent gases in C2 and C3 hydrocarbon products by GC-PDD according to ASTM D8098

This method covers the determination of hydrogen, nitrogen, oxygen, methane, carbon monoxide, and carbon dioxide in the parts per billion mole (nmol/mol) to parts per million mole (µmol/mol) range in C2 and C3 hydrocarbons

#### Instrumentation


The trace gas analyser is configured using a Thermo Trace GC1600 with Pulsed Discharge Detector (PDD or PDHID), 2 purged diaphragm valves and 3 separation columns. The 10 port valve injects the sample and backflushes the matrix compound using column 1.  $CO_2$  is separated on column 2, and elutes to the PDD via restriction R.  $H_2$ ,  $O_2$ ,  $N_2$ ,  $CH_4$  and CO are temporary isolated and separated using column 3 and subsequently detected.

The used columns are:

- 1. Rt-Qbond, 30 m × 0.53 mm
- 2. Rt-Qbond, 30 m × 0.53 mm
- 3. Rt-MSieve 30 m × 0.53 mm

# Very leak-tight diaphragm valves | Low maintenance

Ultra-durable diaphragm valves with an internal purge option offer low maintenance and a very low background level, providing the best possible sensitivity and dynamic range. The valves are mounted in a separately heated valve compartment to ensure stable performance.



**Figure 2** Chromatogram ASTM D8098 analyser




Figure 1 ASTM D8098 Analyser diagram

## **Excellent analytical performance**

The GAS ASTM D8098 analyser shows excellent analytical performance. See the data below and the chromatogram (figure 2).

Repeatability: < 0.5-3% RSD (depending on concentration)


Linearity: > 0.98

Runtime: 7 minutes

Detection limit (s/n=3\*noise):

Н,	< 20 ppb
0,	< 10 ppb
N,	< 10 ppb
CH,	< 10 ppb
co	< 20 ppb
CO <sub>2</sub>	< 10 ppb





**Figure 3** ASTM D8098 analyser is based on Trace GC1600 (left), and is also available on CompactGC<sup>4.0</sup> (right).

## Conclusion

The High Purity Analyser from GAS fully complies with the requirements of ASTM D8098 for the analysis of trace levels of permanent gases in C2 and C3 hydrocarbon products.