

Chromeleon CDS extension

Roundabout

24/7 continuous data acquisition and reporting

AO100 v1.1

Roundabout enhances Chromeleon CDS by automating sample list generation, streamlining analytical workflows, and supporting continuous system operation. It also features integrated auto-calibration and automated summary reporting for improved data accuracy and compliance.

Roundabout software option for infinite data acquisition and reporting

Key features include:

- Automatic generation of sample lists for uninterrupted operation
- Advanced summary reports for comprehensive data reporting and quality charting
- Automatic re-calibration based on system performance
- Result verification against predefined acceptance limits
- Time stamping for each injection
- Tracking of injection counts for audit and maintenance purposes




Figure 1

Roundabout - example of main program window

Ordering information: AO100

Requirements:

- Windows 10 or 11 operating system
- Chromeleon 7.3.2 maintenance update D
- A fully setup and functioning instrument, eWorkflow and methods

About GAS

Global Analyser Solutions provides GC & GC-MS solution for Energy, Refinery, Chemical and Environmental markets. Our analysers address a broad spectrum of measuring requirements with high precision and reliability. Please reach out for more information on our website.